

WATCHMAN™

Left Atrial Appendage Closure Device

Clinical Study Highlights



WATCHMAN™ Clinical – Portfolio Summary

Pilot	➡	Early feasibility with >6 years of follow-up	➡	N=66
PROTECT-AF	➡	Primary efficacy, CV death, and all-cause mortality superior to warfarin at 4 years¹	➡	N=707
CAP Registry	➡	Significantly improved safety results²	➡	N=566
ASAP	➡	Expected rate of stroke reduced by 77% in patients contraindicated to warfarin^{3,4}	➡	N=150
PREVAIL	➡	Improved implant success procedure safety confirmed with new and experienced operators⁵	➡	N=407
CAP2	➡	579 patients enrolled in 48 sites	➡	N=579
Meta-Analysis	➡	Evaluation and analysis of the totality of WATCHMAN trial data⁶	➡	N=1,877
EWOLUTION	➡	Compilation of real-world clinical outcomes data for WATCHMAN^{7,8}	➡	N=1,025

1 Reddy, VY et al. JAMA. 2014; 312(19):1988-1998

2 Reddy, VY et al. Circulation. 2011;123:417-424;

3 Reddy, VY et al. JACC.2013 Jun 25;61(25):2551-6.

4 Sharma D et al. JACC 2016 May 10;67(18):2190-2192.

5 Holmes DR et al. JACC 2014; Jul 8;64(1):1-12

6 Reddy, VY et al. JACC 2017 Dec 19;70(24):2964-2975

7 Boersma LVA et al., Eur Heart J. 2016;37(31):2465-74.

8 Boersma LVA et al., Heart Rhythm. 2017 Sep;14(9):1302-1308.

WATCHMAN™ - Most Studied LAAC Device

Only one proven with long-term data from randomized trials and multi-center registries

	PROTECT AF ¹	CAP Registry ¹	PREVAIL ¹	CAP2 Registry ¹	EWOLUTION ²
Enrollment	2005-2008	2008-2010	2010-2012	2012-2014	2013-2015
Purpose	Demonstrate safety and effectiveness of the WATCHMAN implant compared to long-term warfarin	Continued Access Registry	Demonstrate safety and effectiveness of the WATCHMAN implant compared to long-term warfarin	Continued Access Registry	Collect peri-procedural safety data and long-term outcomes data for patients implanted with WATCHMAN™ in a real world clinical setting
Study Design	2:1 Randomized, non-inferiority	Non-randomized	2:1 Randomized, non-inferiority	Non-randomized	Prospective, single-arm, multi-center, real-life registry
Primary Endpoints	1. Effectiveness: Stroke, systemic embolism and cardiovascular/unexplained death 2. Safety: Life-threatening events, which include device embolization requiring retrieval and bleeding events		1. Effectiveness: Stroke, systemic embolism and cardiovascular/unexplained death 2. Effectiveness: Ischemic stroke or systemic embolism, occurring after 7 days post-randomization or WATCHMAN implant procedure 3. Safety: Death, ischemic stroke, systemic embolism and procedure/device-related complications within 7 days of implantation procedure		Primary analysis includes procedural complications, incidence of stroke, bleeding, and death


1: FDA Oct 2014 Panel Sponsor Presentation.

2: Boersma L.V.A. et al., Catheter Cardiovasc Interv. 2016 Sep;88(3):460-5.

PROCEDURAL SUCCESS AND SAFETY



Implant success results



98.5%
Highest implant success rate
of all WATCHMAN trials¹⁻³


Implant success in EWOLUTION when compared to prior WATCHMAN studies

1. PREVAIL, Holmes, DR et al. JACC 2014;Vol.64, No.1

2. (2014). Circulatory System Devices Panel: WATCHMAN® Left Atrial Appendage Closure Therapy Sponsor Presentation. 2014 FDA Circulatory System Devices Panel.

3. EWOLUTION, Boersma L.V.A et al. Eur Heart J. 2016 Aug;37(31):2465-74


Procedural results & Seal^{1,2}



1: Boersma LV. Et al., EWOLUTION trial: 1-Year Efficacy and Safety of Left Atrial Appendage Closure with WATCHMAN in Patients with or without Contraindication to Oral Anticoagulation.– *Heart Rhythm*. 2017 Sep;14(9):1302-1308

2: Boersma LVA et al., Implant Success and Safety of Left Atrial Appendage Closure with the WATCHMAN device: peri-procedural outcomes from the EWOLUTION Registry. *Eur Heart J*. 2016;37(31):2465-74.

Favorable Procedural Safety Profile



2.8%
Lowest peri-procedural risk
of all WATCHMAN trials¹⁻³

SAE: Serious Adverse Event - Composite of vascular complications includes cardiac perforation, pericardial effusion with tamponade, ischemic stroke, device embolization, and other vascular complications¹

¹Includes observed PE not necessitating intervention, AV fistula, major bleeding requiring transfusion, pseudoaneurysm, hematoma and groin bleeding

1. PREVAIL, Holmes, DR et al. JACC 2014;Vol.64, No.1

2. (2014). Circulatory System Devices Panel: WATCHMAN® Left Atrial Appendage Closure Therapy Sponsor Presentation. 2014 FDA Circulatory System Devices Panel.




3. EWOLUTION, Boersma L.V.A et al. Eur Heart J. 2016 Aug;37(31):2465-74

EFFICACY

Stroke Risk Reduction



PREVAIL Endpoints: Final 5-Year Results


Endpoint	Definition	Final Result
Primary Safety	Acute (7-day) occurrence of death, ischemic stroke, systemic embolism and procedure or device related complications requiring major cardiovascular or endovascular intervention	
Primary Efficacy Composite	Comparison of composite of stroke, systemic embolism, and CV/unexplained death	
Primary Efficacy Ischemic Stroke/SE	Comparison of ischemic stroke or SE occurring >7 days post randomization	

All PREVAIL analyses were pre-specified to use an informative prior that included a portion of PROTECT AF

Source: Reddy VY, et al. 5-Year Outcomes After Left Atrial Appendage Closure: From the PREVAIL and PROTECT AF Trials. JACC 2017 Dec 19;70(24):2964-2975

Patient Level Meta-Analysis

PROTECT AF, PREVAIL 5 Years¹



1: Reddy VY, 5-Year Outcomes After Left Atrial Appendage Closure From the PREVAIL and PROTECT AF Trials. JACC 2017 Dec 19;70(24):2964-2975

Patient Level Meta-Analysis - PROTECT AF, PREVAIL 5 Years¹

Key Outcomes:


At 5 years, WATCHMAN showed:

- **Comparable stroke risk reduction to warfarin**
- **Statistically significant reductions vs. warfarin**




1: Reddy VY, 5-Year Outcomes After Left Atrial Appendage Closure From the PREVAIL and PROTECT AF Trials. JACC 2017 Dec 19;70(24):2964-2975

PROTECT AF Ischemic Strokes: Same Rate Once Accounting For Procedure-related Strokes




Reddy, V. et al. *Circulation* 123(4): 417-424.

Warfarin Ischemic Stroke Rate in PREVAIL Differs from Other Trials



1. Reddy JACC (2017 Dec 19;70(24):2964-2975). 2. Miller. AJC (2012). 3. Giugliano. NEJM (2013)


WATCHMAN Comparable to warfarin for Ischemic Stroke in 2 RCT, 4 Registries



1: Reddy VY, 5-Year Outcomes After Left Atrial Appendage Closure From the PREVAIL and PROTECT AF Trials. JACC 2017 Dec 19;70(24):2964-2975

WATCHMAN Significant Reduction in Disabling Strokes

(Patient-Level Meta-Analysis)



WATCHMAN

warfarin


Disabling Stroke defined as MRS ≥ 2

Two strokes in PREVAIL are excluded because the baseline MRS score was unavailable

1: Reddy VY, 5-Year Outcomes After Left Atrial Appendage Closure From the PREVAIL and PROTECT AF Trials. JACC 2017 Dec 19;70(24):2964-2975

EWOLUTION

Effectiveness in ischemic stroke reduction – annual stroke rate¹



84% reduction in ischemic strokes (annual stroke rate 1.1%) as compared to no therapy*.

*Effectiveness in stroke reduction vs. estimated in the absence of therapy for comparable CHA₂DS₂-VASc scores based on Friberg et al. EHJ 2012

1: Boersma LV. Et al., EWOLUTION trial: 1-Year Efficacy and Safety of Left Atrial Appendage Closure with WATCHMAN in Patients with or without Contraindication to Oral Anticoagulation.– *Heart Rhythm*. 2017 Sep;14(9):1302-1308

EFFICACY

Bleeding Reduction



Bleeding Risks Compound Over Time

CHA ₂ DS ₂ - VASc* Score	Annual % Stroke Risk	HAS-BLED** Score	Annual % Bleed Risk	10-Year Bleeding Risk (%)***
0	0	0	0.9	8.6
1	1.3	1	3.4	29.2
2	2.2	2	4.1	34.2
3	3.2	3	5.8	45.0
4	4.0	4	8.9	60.6
5	6.7	5	9.1	61.5


* 2014 AHA / ACC / HRS Guidelines

** Lip. JACC (2011)

*** Assumes constant risk despite increasing age and bleeding risk is independent from bleeding risk in previous years

EWOLUTION

Effectiveness in bleeding reduction – Annual bleeding rate¹



48% reduction in major bleeding events compared to warfarin.

54% reduction in major bleeding events if excluding periprocedural bleeding.

*Effectiveness in bleeding reduction vs. estimated under VKA therapy for comparable HAS-BLED scores based on Lip et al. JACC 2011

1: Boersma LV. Et al., EWOLUTION trial: 1-Year Efficacy and Safety of Left Atrial Appendage Closure with WATCHMAN in Patients with or without Contraindication to Oral Anticoagulation.– *Heart Rhythm*. 2017 Sep;14(9):1302-1308

72% Major Bleeding Reduction Long Term Post-6 months

Post Procedure Therapy

Destination Therapy

Warfarin +
ASA
(81mg) daily

Clopidogrel (75mg) + ASA
(325 mg) daily

ASA (325mg) daily

Implant

45 days*

6 months


*if leak >5mm, patients remained on warfarin + ASA until seal documented, skipping the clopidogrel + ASA pharmacotherapy

	LAAC (n=732)		Long-term warfarin (n=382)		Rate Ratio	P value
	Bleeding Rate (n events / N at risk)	Event Rate per 100 pt-yrs (n events / N at risk)	Bleeding Rate (n events/N at risk)	Event Rate per 100 pt-yrs (n events / N at risk)		
Overall	10.8 (79/732)	3.5 (79/2268)	11.3 (43/382)	3.6 (43/1187)	0.96	0.84
Post Procedure	5.9 (40/682)	1.8 (40/2255)	11.3 (43/381)	3.6 (43/1180)	0.49	0.001
Destination	3.2 (19/601)	1.0 (19/1958)	9.7 (35/360)	3.5 (35/1004)	0.28	<0.001

Overall period defined as after randomization to the end of follow-up; post-procedural period as >7 days after randomization to the end of follow-up; destination therapy period as beyond 180 days post-randomization, when patients assigned to LAA closure were eligible to receive aspirin alone.

Bleeding Outcomes after Left Atrial Appendage Closure Compared with Long-term Warfarin

Freedom of Major Bleeding Over 3 Adjunctive Pharmacotherapy Intervals



Price, M. J., V. Y. Reddy, et al. JACC: CV Interv 2015; 8(15): 1925-1932

Conclusion

- WATCHMAN LAAC is a proven and safe therapy
- WATCHMAN is the only LAAC device with 5-year efficacy data from 3 major studies
- WATCHMAN demonstrated **comparable stroke risk reduction** and **statistically significant reductions** in **disabling/fatal strokes, major non-procedure related bleeding** and **mortality** compared to warfarin¹:
 - 55% reduction in disabling/fatal stroke
 - 52% reduction in major non-procedure related bleeding
 - 41% reduction in CV/unexplained mortality

1. Reddy VY et al. JACC 2017 Dec 19;70(24):2964-2975